



# NowCheck

Total Solution for COVID-19 Diagnosis



- NowCheck COVID-19 Ag Test
  - Nasopharyngeal swab version
  - Nasal swab version
- NowCheck COVID-19 IgM/IgG Test

# NowCheck COVID-19 Ag Test

The NowCheck COVID-19 Ag Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human nasopharynx or nasal cavity.

## Specifications


- Specimen : Nasopharyngeal swab, Nasal swab
- Testing time : 15~30 min.
- Packing Unit : 25 Tests/Kit
- Storage Condition : 2~30°C (36~86°F)



## Test Procedure


### Nasopharyngeal swab

- 1 Insert a nasopharyngeal swab into the nostril of the patient, and rotate the swab over the posterior nasopharynx surface 3~4 times.
- 2 Insert the swab into an extraction buffer tube, and stir it more than 5 times while squeezing the tube.




### Nasal swab

- 1 After getting a patient to blow one's nose, insert a nasal swab into one of the nostrils, and rotate the swab against the nasal wall 5-6 times. Repeat in the other nostril, using the same swab.
- 2 Insert the swab into an extraction buffer tube, and stir it more than 10 times while squeezing the tube.




### Analysis

- 3 Remove the swab while squeezing the sides of the tube, and press the nozzle cap tightly onto the tube.



- 4 Apply 3 drops of the extracted specimen to the specimen hole of the test device.




- 5 Read the test result after 15~30 min. The test can be read up to 30 min.



\*\* Do not read test results after 30 min. It may give false results.

## Test Result




# NowCheck COVID-19 Ag Test

## Clinical Performance

### :: Evaluation study in Brazil (January – February 2021, FIND)

- A total of 218 nasal & nasopharyngeal swab specimens from symptomatic patients were tested.



Specimens from symptomatic patients (N=218)		RT-PCR (Nasopharyngeal)		
		Positive	Negative	Total
NowCheck COVID-19 Ag (Nasal, Nasopharyngeal)	Positive	71	2	73
	Negative	8	137	145
	<b>Total</b>	<b>79</b>	<b>139</b>	<b>218</b>

	NowCheck COVID-19 Ag (Nasopharyngeal / Nasal)
Clinical Sensitivity (95% CI)	89.9% (81.3, 94.8)
Sensitivity days ≤ 7	92.5% (83.7, 96.8)
Sensitivity Ct ≤ 33	97.2% (90.4, 99.2)
Sensitivity Ct ≤ 25	100% (92.3, 100)
Clinical Specificity (95% CI)	98.6% (94.9, 99.6)

### :: Evaluation study in Brazil (July – August 2020, FIND)

- A total of 400 nasopharyngeal swab specimens from symptomatic patients were tested.

Specimens from symptomatic patients (N=400)		RT-PCR (Nasopharyngeal)		
		Positive	Negative	Total
NowCheck COVID-19 Ag (Nasopharyngeal)	Positive	91	8	99
	Negative	11	290	301
	<b>Total</b>	<b>102</b>	<b>298</b>	<b>400</b>

	NowCheck COVID-19 Ag (Nasopharyngeal)
Clinical Sensitivity (95% CI)	89.2% (81.7, 93.9)
Sensitivity days ≤ 7	92.2% (84.8, 96.2)
Sensitivity Ct ≤ 33	91.4% (83.9, 95.6)
Sensitivity Ct ≤ 25	94.8% (85.9, 98.2)
Clinical Specificity (95% CI)	97.3% (94.8, 98.6)

# NowCheck COVID-19 Ag Test

## SARS-CoV-2 Variant Study

### :: Internal evaluation study (July 2021, BIONOTE Inc. Laboratory)


- *In-silico* analysis of nucleocapsid (N) protein
- Analytical sensitivity tests using recombinant N protein
- Analytical sensitivity tests using cultured virus

Variants		First identified country	<i>In-silico</i> analysis	Analytical sensitivity test	
Name	Lineage			Recombinant protein	Cultured virus
Wuhan-Hu-1	B	China	Confirmed	Confirmed	Confirmed
Alpha	B.1.1.7	United Kingdom		Confirmed	Confirmed
Beta	B.1.351	South Africa		Confirmed	Confirmed
Gamma	P.1	Brazil		Confirmed	To be tested
Kappa	B.1.617.1	India		Confirmed	To be tested
Delta	B.1.617.2	India		Confirmed	To be tested

### :: Evaluation study in Switzerland (June 2021, University of Geneva & FIND)

- Analytical sensitivity tests using cultured virus
- NowCheck COVID-19 Ag Test showed
  - Better (2~16 times) performance on all variants than on the early pandemic lineage
  - Best performance on the Gamma variant among the 9 commercial kits evaluated together

Variants	Early pandemic lineage (B.1.160)	Alpha (B.1.1.7)	Beta (B.1.351)	Zeta (P.2)	Gamma (P.1)
Analytical sensitivity (PFU/mL)	$10^{3.33}$ / Control	$10^{3.03}$ / 2 times	$10^{3.03}$ / 2 times	$10^{2.73}$ / 4 times	$10^{2.13}$ / 16 times



**NowCheck COVID-19 Ag Test provides reliable and highly accurate results for detecting the abovementioned SARS-CoV-2 variants.**

# NowCheck COVID-19 Ag Test

## Analytical Performance

### :: Limit of Detection (LoD)

SARS-CoV-2 strain tested	Virus stock titer	Specimen type	LoD (Spiking titer)	Final working titer
NCCP 43326/2020 /Korea	1 X 10 <sup>6.2</sup> TCID <sub>50</sub> /ml	Direct nasopharyngeal swab	3.12 X 10 <sup>2.2</sup> TCID <sub>50</sub> /ml	6.24 X 10 <sup>1.2</sup> TCID <sub>50</sub> /ml
		Direct nasal swab	1.56 X 10 <sup>2.2</sup> TCID <sub>50</sub> /ml	5.37 X 10 <sup>1.2</sup> TCID <sub>50</sub> /ml

### :: No Cross-reactivity

No cross-reactivity with potential cross-reactive substances except SARS-CoV

Category	Sort	Result	Category	Sort	Result	
Virus	SARS-coronavirus	POS	Bacteria	Haemophilus influenzae	NEG	
	MERS-coronavirus	NEG		Mycoplasma pneumoniae	NEG	
	Human Coronavirus	NEG		Streptococcus pneumoniae	NEG	
	Influenza A	NEG		Streptococcus pyogenes	NEG	
	Influenza B	NEG		Streptococcus salivarius	NEG	
	Respiratory syncytial virus	NEG		Candida albicans	NEG	
	Human Metapneumovirus (hMPV)	NEG		Bordetella pertussis	NEG	
	Parainfluenza virus	NEG		Moraxella catarrhalis	NEG	
	Rhinovirus	NEG		Pseudomonas aeruginosa	NEG	
	Enterovirus	NEG		Staphylococcus epidermidis	NEG	
	Adenovirus	NEG		Mycobacterium tuberculosis	NEG	
	Human immunodeficiency virus lysate	NEG		Others	Pooled human nasal wash – representative of normal respiratory microbial flora	NEG

### :: No Cross-reactivity & Interference

- No interference with endogenous/exogenous interfering substances

### :: Inactivation by the extraction buffer

- SARS-CoV-2 was inactivated by the extraction buffer in 2 minutes.

Type	Virus Spiking	Cytopathic Effect	Interpretation
Extraction buffer	0	No CPE	Virus inactivated
Cell culture media		CPE	Positive control

### :: Matrix Equivalency

- Clinical matrix including swabs and VTM do not affect the detection of COVID-19 Ag.

# NowCheck COVID-19 IgM/IgG Test

NowCheck COVID-19 IgM/IgG Test is a rapid chromatographic immunoassay for the qualitative detection of specific antibodies to SARS-CoV-2 present in human serum, plasma, or whole blood.

## Specifications




- Specimen : Whole blood, Serum, Plasma
- Testing time : 10~15 min.
- Packing Unit : 25 Tests/Kit
- Storage Condition : 2~30°C (36~86°F)



## Test Procedure




### [ Using Capillary whole blood ]

- 1 Collect the 20  $\mu$ l of capillary whole blood to the black line
- 2 Add the collected blood to the specimen hole
- 3 Add 3 drops (90  $\mu$ l) of buffer vertically
- 4 Read the test result in 10~15 min.




### [ Using serum/plasma/whole blood ]

- 1 Collect the serum, plasma (10  $\mu$ l) or whole blood (20  $\mu$ l) with a micropipette
- 2 Add the collected blood to the specimen hole
- 3 Add 3 drops (90  $\mu$ l) of buffer vertically
- 4 Read the test result in 10~15 min.



## Test Result



# NowCheck COVID-19 IgM/IgG Test


## Performance

### :: Clinical evaluation

A total of 105 serum specimens from patients confirmed by Real-time PCR were analyzed. The results are as shown below:


#### - Diagnostic Sensitivity

Positive specimens (N=31)		NowCheck COVID-19 IgM/IgG Test			
		IgM(+)	IgG(+)	IgM(+) and/or IgG(+)	
Days (post-symptom onset)	0 - 7	100% (1/1)	100% (1/1)	100% (1/1)	96.8% (30/31)
	8 - 14	66.7% (2/3)	100% (3/3)	100% (3/3)	
	≥ 15	88.9% (24/27)	96.3% (26/27)	96.3% (26/27)	



#### - Diagnostic Specificity

Negative specimens (N=74)		NowCheck COVID-19 IgM/IgG Test		
		IgM(-)	IgG(-)	IgM(-) and IgG(-)
Confirmed by Real time-PCR		98.6% (73/74)	100% (74/74)	98.6% (73/74)



The NowCheck COVID-19 IgM/IgG Test showed 96.8% of sensitivity and 98.6% of specificity, compared to the Real-time PCR method.

\* Conducted in July 2020 at Department of Laboratory medicine, Seoul National University Bundang Hospital


\*\* Reference: PowerChek™ 2019-nCoV Real-time PCR kit (KCDC approved, FDA EUA)

### :: Limit of Detection (LoD)

The LoD of NowCheck COVID-19 IgM/IgG Test for SN titer was 1:40.


# Total Solution for COVID-19 Diagnosis

## Why Do We Need Antigen and Antibody Tests?



## Diagnostic Algorithm

For individuals suspected to have COVID-19



## Ordering Information

Product No.	Product Name	Product Type	Packing Unit
RG1901DG (Nasopharyngeal) RG1901DGN (Nasal)	NowCheck COVID-19 Ag Test	Device	25 Tests/Kit
RB2901DG	NowCheck COVID-19 IgM/IgG Test	Device	25 Tests/Kit
RG1901CD	NowCheck COVID-19 Ag Control	Control	10 Tests/Kit
RB2901CD	NowCheck COVID-19 IgM/IgG Control	Control	10 Tests/Kit