## ONE STEP Newcastle Disease virus Antigen TEST

For veterinary diagnostic use only

# Anigen Rapid NDV Ag Test Kit

## ■ Principles

The **Rapid NDV Ag Test Kit** is a chromatographic immunoassay for the qualitative detection of Newcastle Disease virus(NDV) antigen in oropharynx, spleen, or kidney from chicken.

The **Rapid NDV Ag Test Kit** shows two letters which are test(T) line and control(C) line on the surface of device. The control line is a reference line which provides the test is performing properly.

The control line has to be appeared at all times after the test has performed. If the NDV antigens exist in sample, a purple test line would appear in the result window.

The highly selective antibodies to NDV are used as a capture and detector in the assay. These antibodies are capable of detecting NDV antigen in chicken oropharynx, spleen or kidney.

## ■ Materials provided

Reagent	10 Tests/Kit
Anigen Rapid NDV Ag Test Device	10
Assay diluent tube	10
Disposable swab	10
Disposable dropper	10
Instructions for use	1

## ■ Materials required, but not provided

- 1) Timer

## ■ Precautions

- 1) The Rapid Test Kit is for chicken only. Do not apply to other animals.
- 2) The test device is sensitive to humidity as well as heat. Perform the test immediately after removing the test device from the foil pouch.
- 3) Do not reuse the test components.
- 4) Apply the assay diluent vertically.
- 5) Do not touch the membrane in the result window of test device.
- 6) Do not use the test kit beyond the stated expiration date marked on the package label.
- 7) Do not use the test kit if the pouch is damaged or the seal is broken.
- 8) Do not mix components from different lot numbers because the components in this kit have been quality control tested as standard batch unit.
- 9) All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- 10) Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.

## ■ Storage and Stability

- 1) Store the test kit at 2~30°C. **DO NOT FREEZE.**
- 2) Do not store the test kit in the direct sunlight.
- 3) The test kit is stable within the expiration date that marked on the package label.


## ■ Collection and Preparation of Sample

- 1) Chicken oropharynx, spleen, or kidney should be used for this test.
- 2) Poke the organ(oropharynx, spleen, or kidney) using a swab. Or Make the 20~50% homogenates of spleen or kidney using PBS. Ex) 50% homogenate: Mix 2g of spleen with 2ml of PBS, 20% homogenate: Mix 1g of spleen with 4ml of PBS.
- 3) If sample is not tested immediately, it should be refrigerated at 2~8°C for up to 24 hours. For longer storage, freeze at or below -20°C in vials with airtight seals.

## ■ Procedure of the Test

- 1) Allow all kit components and samples to room temperature before use.
- 2) Collect the samples using a swab.

- 3) Insert the swab into the assay diluent tube and mix a swab for 10 seconds.

- 4) Wait for 1 minutes to settle down the samples.

- 5) Remove the test device from the foil pouch, and place it on a flat and dry surface.
- 6) Using a disposable dropper provided, take the supernatant sample in the tube.

- 7) **Apply four(4) drops of sample** into the sample hole, drop by drop.

- 8) Start the timer. The sample will flow across the result window. If it does not appear after 1 minute, add one more drop of prepared sample to the sample hole.


- 9) Interpret test results at **10 minutes**. Do not read the result after 20 minutes.

## ■ Interpretation of the Result

### 1) Negative result

One Control line "C" appears in the result window.

### 2) Positive result

Test("T") line and control("C") line within the result window indicate the presence of NDV antigen.

### 3) Invalid Result

If the control("C") line does not appear, the result might be considered invalid. The samples should be retested.

## ■ Limitations of the Test

- 1) Although the Rapid NDV Ag Test kit is very accurate in detecting Newcastle Disease virus antigen, a low incidence of false results can be occurred. Other clinically or laboratory tests might be required if questionable results are obtained. As other diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be diagnosed by the veterinarian after all clinical and laboratory findings have been evaluated. The detection limit of this kit is approximately 0.25 HAU.
- 2) The reading window may show a light pink background coloration; this will not affect the accuracy of the results.
- 3) BioNote and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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