## ONE STEP Bovine *Brucella abortus* Antibody TEST

For veterinary diagnostic use only

### Anigen Rapid B.Brucella Ab Test Kit

#### ■ Principles

The Anigen Rapid B.Brucella Ab Test Kit is a chromatographic immunoassay for the qualitative detection of antibodies against *Brucella abortus* in bovine whole blood, plasma, serum or milk.

The Anigen Rapid B.Brucella Ab Test Kit has two letters which are test(T) line and control(C) line on the surface of device. Test line and control line in the result window are not visible before applying any samples. The control line is a reference line which indicates the test is performing properly. The control line has to appear every time when the test has performed. If the antibodies against *Brucella abortus* are present in sample, a purple test line would appear in the result window.

#### ■ Materials provided

Reagent	10 Tests/Kit
Anigen Rapid B.Brucella Ab Test Device	10
Assay diluent bottle	1
Disposable capillary tube	10
Instructions for use	1

Figure1. A black line on the capillary tube is the indicator line for 10 $\mu$ l.

#### ■ Materials required, but not provided

- 1) Timer

#### ■ Precautions

- 1) The test kit is for bovine use only. Do not use for other animals.
- 2) The test device is sensitive to humidity as well as heat. Perform the test immediately after removing the test device from the foil pouch.
- 3) Do not reuse the test components.
- 4) Apply the sample and assay diluent vertically.
- 5) Do not touch the membrane in the result window of test device.
- 6) Do not use the test kit beyond the stated expiration date marked on the package label.
- 7) Do not use the test kit if the pouch is damaged or the seal is broken.
- 8) Do not mix components from different lot numbers because the components in this kit have been quality control tested as standard batch unit.
- 9) All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- 10) Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.

#### ■ Storage and Stability

- 1) Store the test kit at 2~30°C. **DO NOT FREEZE.**
- 2) Do not store the test kit in the direct sunlight.
- 3) The test kit is stable within the expiration date that marked on the package label.


#### ■ Collection and Preparation of Sample

- 1) Bovine whole blood, serum, plasma or milk should be for this test.  
**[Whole blood]** Collect the whole blood into the anticoagulant tube. If anticoagulated whole blood is not immediately tested, they should be refrigerated at 2~8°C and used within 24 hours.  
**[Serum]** Collect the whole blood into the collection tube (NOT containing anticoagulants such as heparin, EDTA and sodium citrate), leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get supernatant.  
**[Plasma]** Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) and then centrifuge blood to get plasma.  
**[Milk]** Collect the milk into the collection tube.
- 2) If serum, plasma or milk samples are not tested immediately, they should be refrigerated at 2~8°C. For longer storage, freezing

is recommended. Frozen samples should be brought to room temperature (15~30°C) prior to use.

- 3) Samples containing precipitate may yield inconsistent test results. Such samples must be clarified prior to assaying.
- 4) The use of hemolytic, lipaemic, icteric or bacterially contaminated samples should be avoided. Erroneous result may occur.

#### ■ Procedure of the Test

- 1) All reagents and samples must be at room temperature (15~30°C) before use.
- 2) Remove the test kit from the foil pouch, and place it on a flat, dry surface.
- 3) Using a capillary tube, add 10 $\mu$ l of sample into the sample hole (S) on the test device.  

- 4) Add 4 drops (approximately 120 $\mu$ l) of assay diluent into the sample hole(S) vertically.  

- 5) Start the timer. The sample will flow across the result window. If it does not appear after 1 minute, add one more drop of assay diluent.
- 6) Interpret the test results at 20 minutes. Do not interpret after 30 minutes.  

#### ■ Interpretation of the Result


##### 1) Negative

Only control ("C") line appears in the result window.

##### 2) Positive

Test ("T") line and control ("C") line within the result window indicate the presence of antibodies against *Brucella abortus*.

##### 3) Invalid

If the control ("C") line does not appear, the result might be considered invalid. The samples should be retested.

#### ■ Limitations of the Test

- 1) Although the Anigen Rapid B.Brucella Ab Test Kit is very accurate for detecting antibodies against *Brucella abortus*, a low incidence of false results can be occurred. Other clinical and/or laboratory tests might be required if questionable results are obtained. As other diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be diagnosed by the veterinarian after all clinical and laboratory findings have been evaluated.
- 2) The reading window may show a light pink background coloration; this will not affect the accuracy of the results.
- 3) BioNote and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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