ONE STEP Bovine TB Antibody RAPID TEST

For veterinary diagnostic use only

Anigen Rapid Bovine TB Ab Test Kit

Principles

The **Anigen Rapid Bovine TB Test Kit** is a chromatographic immunoassay for the qualitative detection of *Mycobacterium bovis* antibody in bovine whole blood, plasma or serum.

The Anigen Rapid Bovine TB Test Kit shows two letters which are test (T) line and control (C) line on the surface of device. The control line is a reference line which provides the test is performing properly. The control line has to be appeared at all times after the test has performed. If the antibodies against *Mycobacterium bovis* exist in sample, a purple test line would appear in the result window.

The highly selective recombinant *Mycobacterium bovis* antigens are used as a capture and detector in the assay. These antigens are capable of detecting *Mycobacterium bovis* antibodies in sample with high accuracy.

Materials provided

Materials	10 Tests/Kit
Anigen Rapid Bovine TB Ab Test Devices (See figure 1)	10
Disposable capillary tubes (See figure 2)	10
Developing buffer bottle	1
Whole blood diluent bottle	1
Test tubes	10
Anticoagulant tubes	10
Disposable droppers	10
Paper rack for test tubes	1
Instructions for use	1

Figure 1. Test Device

Figure 2. A black line on the capillary tube is the indicator line for 10µl.

Materials required, but not provided

- 1) Timer

Precautions

- 1) The test kit is for cattle use only. Do not use for other animals.
- 2) The test device is sensitive to humidity as well as heat. Perform the test immediately after removing the test device from the foil pouch.
- 3) Do not reuse test components.
- 4) Apply the sample and developing buffer vertically.
- 5) Do not touch the membrane in the result window of test device.
- 6) Do not use the test kit beyond the stated expiration date marked on the package label.
- 7) Do not use the test kit if the pouch is damaged or the seal is broken.
- 8) Do not mix components from different lot numbers because the components in this kit have been quality control tested as standard batch unit.
- 9) All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- 10) Decontaminate and dispose of all samples, reaction kits and potentially contaminated materials safely in accordance with national and local regulations.
- 11) In case of using whole blood specimen, it may cause false negative result compared to serum and plasma specimens because whole blood specimen is diluted. Use of serum or plasma specimen is recommended for more accurate results excluding emergency situation.

Storage and stability

- 1) Store the test kit at 2~30°C. **DO NOT FREEZE.**
- 2) Do not store the test kit in the direct sunlight.
- 3) The test kit is stable within the expiration date that marked on the package label.


Collection and Preparation of Sample

- 1) Whole blood, serum or plasma should be used as a sample for this test.
[Whole blood] Collect the whole blood into the anticoagulant tube (Max.vol.1.5ml) provided. If anticoagulated whole blood is not immediately tested, they should be refrigerated at 2~8°C and used within 24 hours.
[Serum] Collect the whole blood into the collection tube (NOT containing anticoagulants such as heparin, EDTA and sodium citrate), leave to settle for 30 minutes for blood coagulation and then centrifuge to get serum supernatant.
[Plasma] Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) and then centrifuge to get plasma.
- 2) If serum or plasma samples are not tested immediately, they should be refrigerated at 2~8°C. For longer storage, freezing is recommended. Frozen samples should be brought to room temperature (15~30°C) prior to use.
- 3) Samples containing precipitate may yield inconsistent test results. Such samples must be clarified prior to assaying.
- 4) The use of hemolytic or bacterially contaminated samples should be avoided. Erroneous result may occur.

Procedure of the Test


[Serum & Plasma sample]

- 1) All reagents and samples must be at room temperature(15~30°C) before use.
- 2) Remove the test kit from the foil pouch, and place it on a flat, dry surface.
- 3) Using the disposable capillary tube, add 10µl of serum or plasma to the sample hole (marked "S" on the test device) and wait for 1 minute.
- 4) Add 3 drops of the developing buffer in dropping bottle into the developing buffer hole vertically.
- 5) Start the timer. The sample will flow across the result window. If it does not appear after 1 minute, add one more drop of developing buffer to the developing buffer hole.
- 6) Interpret test results at 20 minutes. Do not interpret after 30 minutes.

[Whole blood sample]

- 1) Remove the test kit from the foil pouch, and place it on a flat, dry surface.
- 2) Add 3 drops of the whole blood diluent into the test tube.
- 3) Add 1 drop (30µl) of a whole blood sample using the disposable dropper and mix them for 1 minute.
- 4) Add 10µl of the mixed sample using a capillary tube to the sample hole (marked "S" on the test device) and wait for 1 minute.
- 5) Add 3 drops of the developing buffer into the developing buffer hole.
- 6) Start the timer. The sample will flow across the result window. If it does not appear after 1 minute, add one more drop of developing buffer to the developing buffer hole.
- 7) Interpret test results at 20 minutes. Do not interpret after 30 minutes.

Interpretation of the Test Results


1) Negative result

The only control ("C") line in the result window appears.

2) Positive result

Test ("T") line and control ("C") line within the result window indicate the presence of antibody against Bovine TB.

3) Invalid result

If the control ("C") line does not appear, the result might be considered invalid. The samples should be retested.

Limitation of the Test

- 1) Although this test Kit is very accurate in detecting antibodies against *Mycobacterium bovis* in the specimen, a low incidence of false results can occur. Other clinical and/or laboratory tests might be required if questionable results are obtained. As other diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be diagnosed by veterinarian after all clinical and laboratory findings have been evaluated.
- 2) The reading window may show a light pink background coloration; this will not affect the accuracy of the results.
- 3) BioNote and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

Expected values

Anigen Rapid Bovine TB Test Kit has been compared with a PPD tests. The overall accuracy is greater or equal to 85.5%.

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