One Step Porcine Transmissible Gastroenteritis Virus and Porcine Epidemic Diarrhea Virus Antigen Test

For veterinary diagnostic use only

Anigen Rapid TGE/PED Ag Test Kit

Principles

The Anigen Rapid TGE/PED Ag Test Kit is a chromatographic immunoassay for the qualitative detection of Porcine Transmissible Gastroenteritis (TGE) Virus antigen and Porcine Epidemic Diarrhea Virus antigen in porcine feces.

The Anigen Rapid TGE/PED Ag Test Kit has two letters which are test("T") line and control("C") line on the surface of device. Test line and control line in the result window are not visible before applying any samples. The control line is a reference line which indicates the test is performing properly. The control line has to appear every time when the test has performed. If the TGE and PED virus antigens are present in sample, a purple test line would appear in the result window.

The highly selective antibodies to TGE and PED virus are used as each capture and detector in the assay. These are capable of detecting TGE and PED virus antigen in sample with high accuracy.

Materials provided (10 Tests/Kit)

Reagent	10 Tests/Kit
Anigen Rapid TGE/PED Ag Test Device	10
Assay diluent tube	10
Disposable swab	10
Disposable dropper	10
Instructions for use	1

Materials required, but not provided

- 1) Timer

Precautions


- 1) Test kit is for porcine only. Do not use for other animals.
- 2) The test device is sensitive to humidity as well as heat. Perform the test immediately after removing the test device from the foil pouch.
- 3) Do not re-use the test components.
- 4) Apply the sample using dropper vertically.
- 5) Do not touch the membrane in the result window of test device.
- 6) Do not use the test kit beyond the stated expiration date marked on the package label.
- 7) Do not use the test kit if the pouch is damaged or the seal is broken.
- 8) Do not mix components from different lot numbers because the components in this kit have been quality control tested as standard batch unit.
- 9) All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- 10) Decontaminate and dispose of all samples, reaction kits and potentially contaminated materials safely in accordance with national and local regulations.

Storage and Stability

- 1) Store the test kit at 2~30°C. **DO NOT FREEZE.**
- 2) Do not store the test kit in the direct sunlight.
- 3) The test kit is stable within the expiration date that marked on the package label.


Collection and Preparation of Sample

- 1) The porcine feces should be used for this test.
- 2) The samples should be tested immediately after collection.
- 3) If samples are not tested immediately, they should be refrigerated at 2~8°C. For longer storage, freezing is recommended. Frozen samples should be brought to room temperature (15~30°C) prior to use.
- 4) The amount of fecal swab may affect the results. It is required to follow the swab amount of feces as shown in the picture on the below. Excessive fecal amount may induce a false positive result and slow migration.


Procedure of the Test


- 1) All reagents and samples must be at room temperature (15~30°C) before use.
- 2) Collect porcine feces samples by using a swab.

- 3) Insert the swab into the assay diluent tube. Mix the swab until the sample has been dissolved into the assay diluent (Approximately 10 sec.).

- 4) Wait for 1 minute to settle down the large particles.

- 5) Remove the test device from the foil pouch, and place it on a flat and dry surface.
- 6) Using the disposable dropper, take the supernatant sample in the tube.

- 7) Apply four(4) drops of sample into each sample hole, drop by drop vertically.


- 8) Start the timer. The sample will flow across the result window. If it does not appear after 1 minute, add one more drop of prepared sample to the sample hole.
- 9) Interpret test results at 10 minutes. Do not read the result after 20minutes.

Interpretation of the Result

1) Negative result


Only control("C") line in the result window appears.


2) Positive result

Test("T") line and control("C") line within the result window indicate the presence of TGE and/or PED antigen.


TGE and PED Ag Positive

TGE Ag Positive





PED Ag Positive

3) Invalid Result

If the control("C") line does not appear, the result might be considered invalid. The sample should be re-tested.

Limitations of the Test

- 1) Although the Anigen Rapid TGE/PED Ag Test Kit is very accurate in detecting TGE antigen and PED antigen, a low incidence of false results can be occurred. Other clinically or laboratory tests might be required if questionable results are obtained. As other diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be diagnosed by the veterinarian after all clinical and laboratory findings have been evaluated. The limit of TGE antigen detection of this kit is approximately $3.13 \times 10^{5.0} \text{TCID}_{50}/\text{mL}$. The limit of PED antigen detection is approximately $10^{5.0} \text{TCID}_{50}/\text{mL}$.
- 2) The reading window may show a light pink background coloration; this will not affect the accuracy of the results.
- 3) BioNote and its distributors cannot be held responsible for the consequences of misuse or misinterpretation of the results given by the test.

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