INTENDED USE

The Vcheck Canine Tnl is an *in vitro* diagnostic test kit for the quantitative measurement of cardiac troponin I concentration in canine serum. Troponins are part of the contractile apparatus of cardiac muscle tissue. The troponin complex consists of 3 different proteins (troponin T, C, and I). Of the troponins, cardiac troponin I (Tnl) is a more specific marker for cardiac myocyte injury. Cardiac injury leads to the release of cardiac Tnl into the circulation, where its concentration correlates with the degree of cardiac damage. The BIONOTE Vcheck Canine Tnl is designed to be used only by veterinarians.

PRINCIPLE

The Vcheck Canine Tnl Test Kit is a fluorescent immunoassay for the quantitative measurement of Canine Tnl concentration. The Vcheck Canine Tnl Test Kit uses specific Canine Tnl antibodies that bind to Canine Tnl. The test procedure involves dissolving fluorescent microparticles conjugated to these specific antibodies into assay diluent that reacts with Canine Tnl in the sample. When the sample is applied to the sample hole of the test device, Canine Tnl in the sample migrates along the nitrocellulose membrane and forms complexes with the anti-Canine Tnl antibodies coated on the membrane. As a result, the density of the test line reflects the concentration of Canine Tnl in the sample. The BIONOTE Vcheck Analyzer reads the density of this test line and calculates the Canine Tnl concentration from the calibration curve data. The control line is a reference line that indicates whether the test has been performed correctly.

MATERIALS PROVIDED

Reagent	5 Tests/Kit
① Vcheck Canine Tnl Test device	5
② Assay diluent tube	5
③ Disposable pipette tip	5
④ Instructions for use	1

MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 100 µl pipette

STORAGE AND STABILITY

- Store the test kit at 1~30°C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

PRECAUTIONS

- The test kit is for canine use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the foil pouch.
- Do not reuse test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers, the components in this kit have been quality control tested as a standard batch unit.
- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with

- national and local regulations.
- Discard all used test materials and specimens after processing autoclave for a minimum of 15 minutes at 121°C.
- Hemolyzed samples or serum with precipitate such as fibrin may give erroneous results.
- This product should be used as a reference for screening purposes, and is a reagent that helps in final diagnosis, not for confirmation. The final diagnosis should be made by a professional veterinarian based on the results of this product and other test results and clinical findings.
- Reference ranges may vary by test population, and each laboratory should establish its own reference interval. A veterinarian must conduct a clinical diagnosis by comprehensively considering the results measured using this product, clinical symptoms, and other test results.
- BIONOTE Vcheck analyzer is recommended to use at 15~30°C.

COLLECTION AND PREPARATION OF SAMPLE


- Canine serum should be used with this test.
- [Serum]**
 - Collect the whole blood into a blood collection tube that does not contain anticoagulant. Leave to settle for 30 minutes for blood coagulation and then centrifuge to get a serum supernatant.
 - The samples should be tested immediately after collection. If serum samples are not tested immediately, they can be stored frozen at -20°C or below for up to 2 months. Samples that have been repeatedly frozen and thawed cannot be used.
 - Serum samples stored frozen at -20°C or below should be brought to room temperature (15~30°C) and be completely thawed and mixed via vortexing or pipetting before use.

TEST PROCEDURE

Allow all kit components and sample to reach room temperature (15~30°C) prior to testing.

[Coding]

- Turn on V200 Analyzer and select "Standard Test".

- Remove the test device from the foil pouch. Once the "Insert device" is displayed in the screen, insert the test device.


[Dilution of sample & Measurement of Tnl value]

- Using a 100 µl pipette, draw 100 µl of sample (serum) and add the sample into an assay diluent tube. And then, mix the sample with diluent by pipetting for 5~6 times.


- Add the mixed sample (100 µl) into the sample hole of the test device and press the "Start" button to initiate testing.

* If the time to press "Start" button is delayed, it may affect the test result.

- The V200 analyzer will display the test result on the screen after 10 minutes.
- Remove the test device.

INTERPRETATION OF THE RESULT

- Read the concentration value of Canine Tnl appearing on the display of the BIONOTE Vcheck Analyzer. (0.01 ~ 20 ng/mL)
- If "↓ 0.01 ng/mL" appears on the display, it means the concentration of Canine Tnl in the specimen is less than 0.01 ng/mL.
- If "↑ 20 ng/mL" appears on the display, it means the concentration of Canine Tnl in the specimen is greater than 20 ng/mL.

REFERENCE RANGE

Concentration	Normal Range
< 0.1 ng/mL	Normal Range
0.1 ~ 0.2 ng/mL	Suspected
> 0.2 ng/mL	Abnormal

* A reference value was established by referring to the difference value between the serum of a healthy dog and the serum of a dog with heart disease described in the literature and the company's evaluation. Each laboratory and the veterinarian in charge must comprehensively consider and evaluate the results of tests such as auscultation, chest radiography, and echocardiography as well as the measurement results of this reagent.

SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. Solution: Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. Solution: Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. Solution: Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. Solution: Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. Solution: Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	The measured total hemoglobin is out of the range of 7 to 23 g/dL. Solution: This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely High Total Hemoglobin	The test is invalid. Solution: Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	

Calibration Overdue	The calibration is overdue. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Not Supported Device	Loading a test device that is not supported by the analyzer. Solution: Check whether the test device is manufactured by BioNote, Inc.
EEE	Internal error has occurred. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.

PERFORMANCE CHARACTERISTIC

- Measuring Range**
Canine Tnl concentration can be measured within the range of 0.01~20 ng/mL by using canine serum. To get numerical results for > 20 ng/mL concentration ranges, samples need to be diluted with negative serum and measured again. To calculate the final Canine Tnl concentration, the measured result must be multiplied by the dilution factor.
- Interference**
No interference was observed for each substance up to the concentration presented in the following table.

Interfering substances	Concentration
Hemoglobin	< 150 mg/dL
Intralipid	< 2,500 mg/dL
Cholesterol	< 250 mg/dL
Bilirubin	< 20 mg/dL
Vitamin C	< 100 mg/dL

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