
Vcheck Feline SAA 2.0

FELINE SERUM AMYLOID A

For veterinary use only

INTENDED USE

The Vcheck Feline SAA 2.0 is an *in vitro* diagnostic test kit for the quantitative measurement of Serum Amyloid A (SAA) concentration in feline serum and plasma. In the case of normal individuals, the concentration of SAA is very low in their serum and plasma but it increases, when a pathological phenomenon occurs such as infection, tissue damage or inflammatory disorders.

An increased concentration of SAA is a non-specific reaction to various disease states and can be used to indicate and evaluate immunoinflammatory responses in related diseases, infection, or tissue damage.

The BIONOTE Vcheck Feline SAA 2.0 is designed to be used only by veterinarians.

PRINCIPLE

The Vcheck Feline SAA 2.0 Test kit is a fluorescent immunoassay for the quantitative measurement of Feline SAA concentration.

The Vcheck Feline SAA 2.0 Test kit uses specific anti-Feline SAA antibodies that will bind to Feline SAA. The test procedure involves dissolving fluorescence microparticles conjugated to these specific antibodies into assay diluent that react with Feline SAA in the sample. When the specimen is applied to the sample hole of the test device, Feline SAA in the specimen migrates along the nitrocellulose membrane and forms complexes with the anti-Feline SAA antibodies coated on the membrane. As a result, the density of the test line reflects the concentration of Feline SAA in the serum or plasma. The BIONOTE Vcheck Analyzer reads the density of this test line and calculates the Feline SAA concentration from the calibration curve data. The control line is a reference line which indicates the test has been performed correctly.

MATERIALS PROVIDED

Reagent	5 Tests/Kit	10 Tests/Kit	20 Tests/Kit
① Vcheck Feline SAA 2.0 Test device	5	10	20
② Assay diluent bottle	5	10	20
③ Disposable pipette tip	10	20	40
④ Instructions for use	1	1	1

MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 5 µl pipette
- 100 µl pipette

STORAGE AND STABILITY

- Store the test kit at 1-30°C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

PRECAUTIONS

- The test kit is for feline use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the foil pouch.
- Do not reuse the test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers; the components in this kit have been quality control tested as a standard batch unit.

- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- Decontaminate and dispose of all samples, used kits, and potentially contaminated materials safely in accordance with national and local regulations.

COLLECTION AND PREPARATION OF SAMPLE

- Feline serum or plasma should be used with this test.
- [Serum]** Collect the whole blood using a plain tube **ONLY**. Leave to settle for 30 minutes for blood coagulation and then centrifuge to obtain a serum supernatant.
[Plasma] Collect the whole blood into a blood collection tube containing anticoagulant (**ONLY** heparin). Then centrifuge to obtain plasma supernatant.
- If serum or plasma samples are not tested immediately, they should be refrigerated at 2-8°C and used within 2 weeks. For longer storage, serum or plasma can be frozen (-20°C or colder). Frozen samples should be brought to room temperature (15-30°C) prior to use.

TEST PROCEDURE

All reagents and samples must be at room temperature (15-30°C) before use.

[Coding]

- Turn on V200 Analyzer and select "Standard Test".


- Remove the test device from the foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.
- After checking Feline SAA item name and test procedure on the display window, proceed as follows.

[Dilution of sample & Measurement]


- Using a 5 µl pipette, draw 5 µl of sample (serum or plasma) and add the sample into an assay diluent bottle (2mL).


- Close the bottle cap and shake for 5-6 times to mix the sample and the dilution thoroughly.

- Add the mixed sample (100 µl) into the sample hole of the test device using the 100 µl pipette and press the [START] to initiate testing.

- Read the concentration value appearing on the display after 5 minutes.
- Remove the test device. The display then shows that BIONOTE Vcheck Analyzer is ready to test the next device.

INTERPRETATION OF THE RESULT

- Read the concentration value of Feline SAA appearing on the display of the BIONOTE Vcheck Analyzer.
- If "Low" appears on the display, it means the concentration of Feline SAA in the specimen is less than or equal to 5 µg/ml.
- If "High" appears on the display, it means the concentration of Feline SAA in the specimen is greater than or equal to 200 µg/ml.

REFERENCE RANGE

< 5 µg/ml	5~10 µg/ml	> 10 µg/ml
Normal	Equivocal	Abnormal

SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. Solution: Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. Solution: Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. Solution: Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. Solution: Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. Solution: Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	The measured total hemoglobin is out of the range of 7 to 23 g/dL. Solution: This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely High Total Hemoglobin	The test is invalid. Solution: Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The calibration is overdue. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	A test device that is not supported by the analyzer has been loaded. Solution: Check whether the test device is manufactured by BioNote, Inc.
Not Supported Device	Internal error has occurred. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
EEE	

PERFORMANCE CHARACTERISTIC

- Measuring Range**
5~200 µg/ml of Feline SAA is measurable in serum and plasma. For > 200 µg/ml results, samples may be diluted with negative serum and re-measured for a numerical result. To calculate the final Feline SAA concentration, the measured result must be multiplied by the dilution factor.
- Interference**
Measurement was not influenced by interfering substances such as triglyceride, bilirubin and anticoagulants (only heparin).
- Analytical specificity**
No other feline blood components are known to cross react with SAA in the Vcheck Feline SAA 2.0 Test kit.

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