

# Comparative evaluation of Vcheck M Ehrlichia/Anaplasma with Real-time PCR

**Key Words** : Vcheck M, Ehrlichia, Anaplasma, Antibody rapid test, Real-time PCR, Sequencing

## Introduction

*Ehrlichia* spp. and *Anaplasma* spp. are rickettsia bacteria that are transmitted by ticks to dogs and cats and are the causative agents of Ehrlichiosis and Anaplasmosis, respectively. They are important tick-borne diseases with a worldwide distribution. They can cause a variety of signs from none to fever and generalized achiness to possible fatality.

The veterinarian usually uses serologic tests and PCR tests to diagnose Ehrlichiosis & Anaplasmosis. In the past, PCR test samples had to be sent to an outside laboratory, but with the Vcheck M, PCR testing is possible in the veterinary clinic.

## Purpose

The goal of this study is to evaluate the diagnostic sensitivity and specificity of the newly developed Vcheck M Ehrlichia/Anaplasma (POCT PCR kit) to laboratory-based real-time PCR.

Additionally, if there are discrepancies between the Vcheck M and lab-based real-time PCR results, sequencing will be requested to an external laboratory to confirm the final result.

## Materials and Methods

SD Biosensor Inc., MDx R&D Department performed tests with Vcheck M and "P" kit (UK) with real-time PCR. Additionally, the sequencing tests were requested to the 'B' Laboratory (Korea) for confirmation.

The samples were a total of 35 canine whole blood checked positive by an antibody rapid test at a laboratory (Paraguay) and an animal hospital (Malaysia).

## Results

The test results for the comparison of Vcheck M and lab-based real-time PCR are described in Tables 1, 2, and 3.

## Conclusion

In this study, there were several discrepancies between the antibody rapid test and the two PCR tests. It is assumed that the dogs were infected with *Ehrlichia* or *Anaplasma* and then recovered.

In addition, discrepancies between Vcheck M and lab-based real-time PCR occurred in a total of 4 samples (Vcheck M positive, PCR negative). As a result of sequencing, these samples were confirmed to be positive.

Based on the results, it was confirmed that Vcheck M Ehrlichia/Anaplasma is superior to existing lab-based PCR in terms of not only convenience but also clinical performance.

<i>Ehrlichia</i> spp.		Reference method (Real-time PCR & sequencing)		
		Positive	Negative	Total
Vcheck M Ehrlichia/Anaplasma	Positive	18	0	18
	Negative	0	17	17
	Total	18	17	35
	Sensitivity	100% (18/18)		
	Specificity	100% (17/17)		

<i>Ehrlichia</i> spp.		Reference method (Real-time PCR & sequencing)		
		Positive	Negative	Total
Lab-based real-time PCR	Positive	15	0	15
	Negative	3	17	20
	Total	18	17	35
	Sensitivity	83.8% (15/18)		
	Specificity	100% (17/17)		

**Table. 1** Sensitivity and specificity of Vcheck M Ehrlichia/Anaplasma and existing lab-based real-time PCR for *Ehrlichia* spp.

<i>Anaplasma</i> spp.		Reference method (Real-time PCR & sequencing)		
		Positive	Negative	Total
Vcheck M Ehrlichia/Anaplasma	Positive	3	0	3
	Negative	0	32	32
	Total	3	32	35
	Sensitivity	100% (3/3)		
	Specificity	100% (32/32)		

<i>Anaplasma</i> spp.		Reference method (Real-time PCR & sequencing)		
		Positive	Negative	Total
Lab-based real-time PCR	Positive	1	0	1
	Negative	2	32	34
	Total	3	32	35
	Sensitivity	33.3% (1/3)		
	Specificity	100% (32/32)		

**Table. 2** Sensitivity and specificity of Vcheck M Ehrlichia/Anaplasma and existing lab-based real-time PCR for *Anaplasma* spp.

Sample No.	<i>Ehrlichia</i> spp.				<i>Anaplasma</i> spp.			
	Antibody Rapid	Vcheck M	Real-time PCR	Sequencing	Antibody Rapid	Vcheck M	Real-time PCR	Sequencing
22	<i>Ehrlichia</i> (+)	<i>Ehrlichia</i> (+)	(-)	<i>E. ewingii</i> (+)	(-)	(-)	(-)	(-)
24	<i>Ehrlichia</i> (+)	<i>Ehrlichia</i> (+)	(-)	<i>E. ewingii</i> (+)	(-)	<i>Anaplasma</i> (+)	(-)	<i>A. platys</i> (+)
27	<i>Ehrlichia</i> (+)	(-)	(-)	(-)	(-)	<i>Anaplasma</i> (+)	(-)	<i>A. platys</i> (+)
29	<i>Ehrlichia</i> (+)	<i>Ehrlichia</i> (+)	(-)	<i>E. canis</i> (+)	(-)	(-)	(-)	(-)

**Table. 3** Sequencing results of discrepant samples between Vcheck M and existing lab-based real-time PCR