Vcheck

D-dimer

CANINE D-DIMER

For veterinary use only

INTENDED USE

The BIONOTE Vcheck D-dimer is an in vitro diagnostic test kit for the quantitative measurement of D-dimer concentration in canine plasma. Since the kit provides quantitative measurement of canine plasma D-dimer levels, BIONOTE Vcheck D-dimer can be used with every dog as a diagnostic tool for canine thromboembolic disorders. D-dimer plasma concentrations can be used to rapidly detect the thrombotic complications and DIC associated with many systemic diseases. The BIONOTE Vcheck D-dimer is designed to be used only by veterinarians.

PRINCIPLE

The Vcheck D-dimer Test kit is a fluorescent immunoassay for the quantitative measurement of canine D-dimer concentration. The Vcheck D-dimer Test kit uses specific Canine D-dimer antibodies that bind to Canine D-dimer. The test procedure involves dissolving fluorescence micro particles conjugated to these specific antibodies into assay diluent that react with Canine D-dimer in the sample. When the sample is applied to the sample hole of the test device, Canine D-dimer in the sample migrates along the nitrocellulose membrane and forms complexes with the anti-D-dimer antibodies coated on the membrane. As a result, the density of the test line reflects the concentration of Canine D-dimer in the sample. The BIONOTE Vcheck Analyzer reads the density of this test line and calculates the Canine D-dimer concentration from the calibration curve data. The control line is a reference line which indicates the test has been performed correctly.

MATERIALS PROVIDED

Reagent	10 Tests/Kit
① Vcheck D-dimer Test device	10
② Assay diluent tube	10
③ Disposable pipette tip	30
④ 1.5 mL tube	10
⑤ Sodium Citrate(1 mL/bottle)	1
⑥ Instructions for use	1

MATERIALS REQUIRED, BUT NOT PROVIDED

1. BIONOTE Vcheck Analyzer
2. 5 µℓ pipette
3. 50 µℓ pipette
4. 100 µℓ pipette
5. A centrifuge

STORAGE AND STABILITY

1. Store the test kit at 2~8 °C. **DO NOT FREEZE.**
2. Do not store the test kit in the direct sunlight.
3. The test kit is stable until the expiry date that is marked on the package label.

Reagent	Open status	Storage	Stability	Note
Test device	Unopened	2~8 °C, Sealed	12 months	Finished product
	Opened	Do not store	-	Use directly
Assay diluent	Unopened	2~8 °C, Sealed	12 months	Finished product
	Opened	Do not store	-	Use directly

PRECAUTIONS

1. The test kit is for canine use only. Do not use for other animals.
2. This reagent needs to be stored at 2~8 °C. If refrigerated, allow all kit components to reach room temperature prior to testing.
3. The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the aluminum foil pouch.
4. Do not reuse test components.
5. Do not touch the membrane in the result window of the test device.
6. Do not use the test kit beyond the stated expiry date marked on the label.
7. Do not use the test kit if the pouch is damaged or the seal is broken.
8. Do not mix components from different lot numbers, the components in this kit have been quality control tested as a standard batch unit.

9. All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
10. Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.
11. Strictly follow the test procedure (E.g. adequate sample volume), as failure to do so may adversely affect test performance and/or produce invalid results.
12. This reagent is designed for quantification of the D-dimer concentration in canine blood using a simple and quick method, but there may be differences in accuracy compared to other laboratory methods.
13. Severely hemolyzed samples or serum with precipitate such as fibrin may give erroneous results.
14. Final diagnosis must be confirmed by a veterinarian with other clinical data available.
15. Use plain tubes only. Failure to do so may adversely affect test performance and/or produce invalid results.
16. BIONOTE Vcheck analyzer is recommended to use at 15~30 °C.

COLLECTION AND PREPARATION OF SAMPLE

1. Canine plasma should be used with this test.
2. [Plasma]

① Using a 50 µℓ pipette, take 50 µℓ of Sodium Citrate and add it into a 1.5 mL tube.

50µℓ

Sodium Citrate

② Add whole blood 450 µℓ into a 1.5 mL tube containing Sodium Citrate (Fill up to 0.5 mL line of 1.5 mL tube)

0.5

0.1

0.5 mL

③ Mix the tube gently using a wrist snap making an 8-character shape. And centrifuge at 6,000 RPM for 5 minutes to separate plasma.

④ If plasma samples are not tested immediately, they should be refrigerated at 2~8 °C and used within 2 weeks. For longer storage, plasma can be frozen (-20 °C or colder). Frozen samples should be brought to room temperature (15~30 °C) prior to use.

\* 1) After collecting blood, immediately separate plasma by mixing Sodium Citrate (3.2% recommended) and whole blood with 1:9 ratio and centrifuge. (At this time, the measurement value might be changed according to the ratio of anticoagulants and whole blood.)

\* 2) If the fresh blood doesn't mix up well and quickly, partial coagulation may occur. Also, decrease of coagulation factors may lead to errors in test results.

TEST PROCEDURE

- \* Allow all kit components and samples to reach room temperature (15~30 °C) for at least 30 minutes prior to testing.
- \* Prepare the same number of the test device, anticoagulants, assay diluent and 1.5 mL tube as the number of the samples.

[Coding]

① Turn on V200 Analyzer and select “Standard Test”.

- ② Remove the test device from the aluminum foil pouch. Once the “Insert Device” is displayed in the screen, insert the test device.

Before starting the test, follow the protocol for preparation of blood plasma (Collection and Preparation of Sample '①', '②', '③', '④'.)

[Dilution of sample & Measurement]

① Using a 5 µℓ pipette, draw 5 µℓ of plasma and add it into an assay diluent tube.

5µℓ

② Use a 100 µℓ pipette to mix the sample with diluent by pipetting for 5~6 times.

100µℓ

③ Add 100 µℓ of the mixed sample into the sample hole of the test device and press the [START] to initiate testing.

100µℓ

④ The V200 analyzer will display the test result on the screen after 5 minutes.

05:00

⑤ Remove the test device.

\* Strictly follow the test procedure including the amount of sample (5 µℓ) used and the test time (5 min), as failure to do so may adversely affect test performance and/or produce invalid results.

INTERPRETATION OF THE RESULT

1. Read the concentration value of Canine D-dimer appearing on the display of the BIONOTE Vcheck Analyzer. (0.1~10 µg/mL)
2. If "↓ 0.1 µg/mL" appears on the display, it means the concentration of Canine D-dimer in the specimen is less than 0.1 µg/mL.
3. If "↑ 10 µg/mL" appears on the display, it means the concentration of Canine D-dimer in the specimen is greater than 10 µg/mL.
4. If the [Invalid] result appears on the screen, a retest shall be carried out.

REFERENCE RANGE

< 0.3 µg/mL	≥ 0.3 µg/mL
Normal	Abnormal TED/DIC* probable

\* TED: Thromboembolic disease,  
DIC: Disseminated intravascular coagulation

PERFORMANCE CHARACTERISTIC

1. Measuring Range  
0.1~10 µg/mL of Canine D-dimer is measurable in plasma (Sodium Citrate). For ↑10 µg/mL results, the sample may be diluted with negative sample and re-measured for a numerical result. To calculate the final Canine D-dimer concentration, the measured result must be multiplied by the dilution factor.
2. Interfering substances study  
No interference was observed for each substance up to the concentration presented in the following table.

Interfering substances	Concentration
Hemoglobin	300 mg/dL
Cholesterol	400 mg/dL
Intralipid	1,000 mg/dL
Bilirubin (total)	15 mg/dL
Vitamin C	62.5 mg/dL

SCREEN MESSAGES AND TROUBLE SHOOTING

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. <b>Solution:</b> Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. <b>Solution:</b> Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. <b>Solution:</b> Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. <b>Solution:</b> Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed.
Barcode Error	<b>Solution:</b> Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely High Total Hemoglobin	The measured total hemoglobin is out of the range of 7 to 23 g/dL. <b>Solution:</b> This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The test is invalid. <b>Solution:</b> Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	The calibration is overdue. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Not Supported Device	A test device that is not supported by the analyzer has been loaded. <b>Solution:</b> Check whether the test device is manufactured by BioNote, Inc.
EEE	Internal error has occurred. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.

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