

FeliCheck-3

General Information

Purpose : Detection of Feline Immunodeficiency virus antibody, Feline Leukemia Virus antigen
Feline Heartworm antibody

Principle : Immunochromatographic assay

Specimen : Whole blood, plasma or serum

Materials provided :


- Anigen Rapid FeliCheck-3 test device
- Assay diluents bottle
- Anticoagulant tube
- Disposable capillary tube

Special Features :

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RC2210CC	FeliCheck-3	Triple Device	5 Tests/Kit