

Anigen Rapid

CDV/CAV Ag

General Information

Purpose : Detection of Canine Distemper virus and Adenovirus antigen

Principle : Immunochromatographic assay

Specimen : CDV Ag - Conjunctiva, CAV Ag - Nasal fluid

Materials provided :


- Anigen Rapid CDV/CAV test device
- Assay diluents tube
- Disposable swab
- Disposable dropper

Special Features :

- Simultaneous detection of CDV and CAV with one specimen preparation
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RC1108DD	CDV/CAV Ag	Dual Device	10 Tests/Kit