

# Caniv-4 (Leish)

## General Information

**Purpose :** Detection of *Dirofilaria immitis* antigen, *Ehrlichia canis* antibody, *Leishmania infantum* antibody *Anaplasma phagocytophilum*/*Anaplasma platys* antibody

**Principle :** Immunochromatographic assay

**Specimen :** Whole blood, plasma or serum

**Materials provided :**


- Anigen Rapid Caniv-4(Leish) test device
- Assay diluents bottle
- Anticoagulant tube
- Disposable capillary tube

**Special Features :**

- Simultaneous detection of four kinds of Vector Borne diseases
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



## Test Procedure



## Ordering Information

Cat. No.	Description	Type	Pack size
RC2127DD	Caniv-4 (Leish)	Quattro Device	10 Tests/Kit