

# CIRD-3 Ag

## General Information

**Purpose :** Detection of Canine Infectious Respiratory Disease (CDV, CAV and CIV)

**Principle :** Immunochromatographic assay

**Specimen :** CDV Ag - Conjunctiva, CAV & CIV Ag - Nasal fluid

**Materials provided :**



- Anigen Rapid CIRD-3 test device
- Assay diluents tube
- Disposable swab
- Disposable dropper

**Special Features :**


- Simultaneous detection of CDV, CAV and CIV with one specimen preparation
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



## Test Procedure



**Interpretation**



## Ordering Information

Cat. No.	Description	Type	Pack size
RC1109DD	CIRD-3 Ag	Dual Device	10 Tests/Kit