

FCoV Ag/Ab

General Information

Purpose : Detection of Feline Coronavirus antigen / antibody

Principle : Immunochromatographic assay

Specimen : Ag - Feline feces
Ab - Whole blood, plasma or serum

Materials provided :

- Anigen Rapid FCov Ag/Ab test device
- Sample collection tube containing assay diluents
- Disposable swab
- Disposable dropper
- Assay diluent dropping bottle
- Anticoagulant tube
- Disposable capillary tube




Special Features :


- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity

Test Procedure

FCoV Ag Test



FCoV Ab Test



Ordering Information

Cat. No.	Description	Type	Pack size
RC1207CC	FCoV Ag/Ab	Dual Device	5 Tests/Kit