

Toxoplasma Ab for Farm

Toxoplasma gondii antibody

Toxoplasmosis is a protozoa disease caused by *Toxoplasma gondii*. The pathogen infects most genera of warm-blooded animals, including humans, but the primary host is the cat family. Infection occurs by eating infected meat, particularly swine products or ingesting water, soil, or food that has come into contact with infected animals' fecal matter. The Anigen Rapid Toxoplasma Ab Test Kit is a chromatographic immunoassay for the qualitative detection of antibody against *Toxoplasma gondii* in serum, plasma or whole blood.




Indications

- Detection of *Toxoplasma gondii* antibody
- Field monitoring of *Toxoplasma* infection
- Sheep and Livestock animal surveillance

Special Features

- One-step testing procedure
- No cross-reaction with other microorganism sera.
- Applicable to various species
- Quick test in 10 minutes
- High sensitivity and specificity
- Sensitivity : 99.9% vs. ELISA
- Specificity : 98.7% vs. ELISA

Test Procedures



Ordering Information

| Cat. No. | Description | Type | Packing size |
|----------|------------------------------|--------|-----------------|
| RB2805DD | Rapid Toxoplasma Ab for Farm | Device | 1 Test x 10/Kit |