

FCoV Ag

General Information

Purpose : Detection of Feline Coronavirus antigen

Principle : Immunochromatographic assay

Specimen : Feline feces

Materials provided :

- Anigen Rapid FCoV test device
- Assay diluents tube
- Disposable swab
- Disposable dropper

Special Features :

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
 - Sensitivity : 92.8 % vs. RT-PCR
 - Specificity : 98.1 % vs. RT-PCR







▪ **Feline infectious peritonitis ; FIP**

There are two basic forms of FIP-effusive (wet) and non effusive (dry) that have been characterized. Cats with effusive FIP usually only survive a few days-weeks at best. Cats with noneffusive FIP can survive for weeks or months, although after neurologic signs begin death usually occurs rapidly.

Test Procedure

[Proper sample amount collected by a swab]



Ordering Information

Cat. No.	Description	Type	Pack size
RG1205DD	FCoV Ag	Device	10 Tests/Kit