

Total IgE (Allergy indicator test)

General Information

Purpose : Detection of Canine Total IgE

Principle : Immunochromatographic assay

Specimen : serum

Materials provided :


- Anigen Rapid Total IgE test device
- Assay diluents bottle
- Disposable capillary tube
- Disposable dropper

Special Features :

- Canine allergy indicator test
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
 - Sensitivity : 97.1% vs. Dog IgE ELISA Quantitation set(BETHYL Lab., Inc.)
 - Specificity : 96.4% vs. Dog IgE ELISA Quantitation set(BETHYL Lab., Inc.)
- Detects serum total IgE higher than 10³U/ml



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RB2110DD	Total Ig E	Device	10 Tests/Kit