

# Caniv-4

## General Information

**Purpose :** Detection of *Dirofilaria immitis* antigen, *Ehrlichia canis* antibody, *Borrelia burgdoferi* antibody *Anaplasma phagocytophilum*/*Anaplasma platys* antibody

**Principle :** Immunochromatographic assay

**Specimen :** Whole blood, plasma or serum

**Materials provided :**


- Anigen Rapid Caniv-4 test device
- Assay diluents bottle
- Anticoagulant tube
- Disposable capillary tube



**Special Features :**

- Simultaneous detection of four kinds of Vector Borne diseases
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity

## Test Procedure



## Ordering Information

Cat. No.	Description	Type	Pack size
RC2120DD	Caniv-4	Quattro Device	10 Tests/Kit